

# Package-Ref <sup>TM</sup> MRD Cocktail Reference Standard CBP90040

### BACKGROUND

Minimal residual disease or Measurable Residual **Disease or Molecular Residual Disease (MRD)** refers to the small number of cancer cells that remain in the body after treatment. The most widely used tests are flow cytometry, polymerase chain reaction (PCR) and nextgeneration sequencing (NGS).

There are two main NGS technical routes for MRD detection of ctDNA: *tumor-agnostic -assays* and tumor-informed assays. Tumor-agnostic -assays only detect plasma, using fixed panels (usually driver genes and targeted drug genes, and multi-omics methods) and analytical methods to detect and analyze ctDNA. Tumorinformed assays require sequencing of primary tumor tissue (usually WES) to identify the patient's specific genomic variation map, and then customize a personalized panel to test ctDNA, which means that both tissue and plasma need to be tested.

#### INTRODUCTION

For ctDNA MRD diagnostic performance evaluation, cb-gene has developed an ultra-low frequency ctDNA standard by Ultrasonic interruption, It covers a wide range of genes and sites and has rich mutation types.

Catalog	Name	Quantity	ſ
CBP90040-1	Package-RefTM MRD Cocktail Reference Standard-0.5%	500 ng/vial	Optimi ddPCR
CBP90040-2	Package-RefTM MRD Cocktail Reference Standard-0.05%	500 ng/vial	Dilutio confirm
CBP90040-3	Package-RefTM MRD Cocktail Reference Standard-0.005%	500 ng/vial	Dilutio confirm
CBP90040-4	Package-RefTM MRD Cocktail Reference Standard-0%	500 ng/vial	ddPCR

**Serial Dilution** 0.005% 0.0005%

GENE

EGFR

EGFR

TP53

CDKN2A

HLA-A

#### How to verify %AF which is ultra-low?

CB-gene adopts methods 1 and 4.

- 2. ultra-Deep sequencing of NGS;
- 3. To check the dilution-fold with the non-human sequence; 4. To check the dilution-fold with the high-%AF mutations.

ized assay on foldnation on foldnation

assay

	EGFR	p.L861Q	c.2582T>A
	ERBB2	p.A775_G776insYVMA	c.2324_2325ins12A7
	KRAS	p.G12V	c.35G>T
	PIK3CA	p.E542K	c.1624G>A
	KRAS	p.Q61H	c.183A>C
	EGFR	p.G719C	c.2155G>T
	ROS1	CD74-ROS1	chr5:149783684:- (h chr6:117646732:- (h
	NRAS	p.Q61K	c.181C>A
	IDH1	p.R132C	c.394C>T
	FBXW7	p.R505C	c.1513C>T
	KMT2C	p.F4496Lfs*21	c.13488delA
	RAC1	p.N92I	c.275A>T
	BCOR	p.Q1321*	c.3961C>T
	BRCA2	p.M784V	c.2350A>G
	ARID1A	p.G1848Wfs*6	c.5540dupG
	CDH1	p.P126Rfs*89	c.377delC
	BARD1	p.K208Rfs*4	c.623delT
	BAX	p.E41Rfs*19	c.121delG
MRD-0.5%	NBN	p.R466Gfs*18	c.1396delT
ctDNA	KMT2D	p.Q809Rfs*121	c.2425delG
	KMT2D	p.P647Hfs*283	c.1940delG
	QKI	p.K134Rfs*14	c.401delA
	KMT2B	p.W336*	c.1007G>A
	NSD1	p.M1531Cfs*43	c.4591delA
	AXIN1	p.E31*	c.91G>T
	AXIN1	p.G508Vfs*197	c.1523del
	ARID1B	p.A547Sfs*35	c.1638 1647delGGC
	HLA-C	p.V100Efs*38	c.299 300delTG
	TP53BP1	p.R1447Efs*36	c.4339delC
	POLD1	N/A	c.E24-2AG>A
	CYLD	p.N722Mfs*13	c.2165delA
	INHA	p.G66Afs*61	c.197delG
	RASA1	p.N838Mfs*4	c.2513delA
	RECOL4	p.G952Afs*92	c.2855delG
	BRAF	p.V600E	c.1799T>A
	ACVR2A	p.K437Rfs*5	c.1310delA
	PIK3CA	p.H1047R	c.3140A>G
	BRCA1	p.D435Y	c.1303G>T
	BRCA2	p.N1784Tfs*7	c.5351delA
	BRCA2	p.E2292A	c.6875A>C

N/A

p.R80\*

p.R45Afs\*32

Mut.

p.S768I

p.E746 A750delELREA

Dilution flod verification	GENE	Mutation	AF% of MRD-0.5% ctDNA	AF% of MRD-0.05% ctDNA	AF% of MRD-0.005% ctDNA
	BRAF	p.V600E	5.33	0.54	0.047
	PIK3CA	p.H1047R	2.81	0.26	N/A
	od BRCA1	p.D435Y	2.53	0.24	N/A
	BRCA2	p.N1784Tfs*7	2.35	0.23	N/A
	TP53	c.783-2A>C	1.02	0.11	N/A
	CDKN2A	p.R80*	0.92	0.1	N/A
	HLA-A	p.R45Afs*32	1.25	0.11	N/A

## **PRODUCT DATA**

- 1. Improve the performance of ddPCR assay;

CDS change	Variant Classification	AF%
c.2236_2250del15GAATTAAGAGAAGCA>-	Pathogenic	0.69
c.2303G>T	Pathogenic	0.61
c.2582T>A	Pathogenic	0.61
c.2324 2325ins12ATACGTGATGGC	Pathogenic	0.41
c.35G>T	Pathogenic	0.39
c.1624G>A	Pathogenic	0.62
c.183A>C	Pathogenic	0.49
c.2155G>T	Pathogenic	0.5
chr5:149783684:- (hg19)		0.55
chr6:117646732:- (hg19)	Pathogenic	0.55
c.181C>A	Pathogenic	0.41
c.394C>T	Pathogenic	0.43
c.1513C>T	Pathogenic	0.75
c.13488delA	Likely Pathogenic	0.46
c.275A>T	Pathogenic	0.4
c.3961C>T	Likely Pathogenic	0.68
c.2350A>G	Benign	0.64
c.5540dupG	Likely Pathogenic	0.61
c.377delC	Pathogenic	0.74
c.623delT	Likely Pathogenic	0.38
c.121delG	Pathogenic	0.75
c.1396delT	Pathogenic	0.71
c.2425delG	Likely Pathogenic	0.51
c.1940delG	Pathogenic	0.44
c.401delA	Likely Pathogenic	0.63
c.1007G>A	Likely Pathogenic	0.75
c.4591delA	Likely Pathogenic	0.74
c.91G>T	Likely Pathogenic	0.66
c.1523del	Likely Pathogenic	0.69
c.1638_1647delGGGCGGCGGC	Pathogenic	0.47
c.299_300de1TG	Likely Pathogenic	0.51
c.4339delC	Likely Pathogenic	0.57
c.E24-2AG>A	Likely Pathogenic	0.73
c.2165delA	Likely Pathogenic	0.57
c.197delG	Likely Pathogenic	0.74
c.2513delA	Likely Pathogenic	0.51
c.2855delG	Likely Pathogenic	0.54
c.1799T>A	Pathogenic	5.33
c.1310delA	Likely Pathogenic	3.98
c.3140A>G	Pathogenic	2.81
c.1303G>T	Uncertain significance	2.53
c.5351delA	Pathogenic	2.35
c.6875A>C	Uncertain significance	1.35
c.783-2A>C	Pathogenic	1.02
c.238C>T	Pathogenic	0.92
c.133del	Likely Pathogenic	1.25

The ddPCR detection system accurately calibrated the AF%.

ACVR2A	BRCA1	HLA-A	KRAS	RASA1
ARID1A	BRCA2	HLA-B	NBN	RECQL4
ARID1B	CDH1	HLA-C	NRAS	ROS1
AXIN1	CDKN2A	IDH1	NSD1	TP53
BARD1	CYLD	INHA	PIK3CA	TP53BP1
BAX	EGFR	KMT2B	POLD1	
BCOR	ERBB2	KMT2C	QKI	
BRAF	FBXW7	KMT2D	RAC1	

30+ genes 40+ loci (black: 0.5% verified genes, red: dilution multiple verification genes)

#### **NOTES:**

- interference from LOB;
- after multiple tests);
- 6,000;

Email:sales@cobioer.com

selected 8 genes and 9 loci as the dilution fold verification of 0.05% selected 2 genes and 2 loci as the dilution fold verification of 0.005%.

For AF=0.5%, each site was accurately calibrated using optimized ddPCR, with an average single-well copy number greater than 6000 to ensure true positivity and distinguish

II. For AF=0.05%, the high-frequency sites (e.g. 5% of BRAF) p.V600E) were used, and the frequency changes of the highfrequency sites were tested by ddPCR (e.g. BRAF p.V600E) diluted 10 fold is 0.5%) to determine the accuracy of the dilution fold, thereby *indirectly verifying* the frequency changes of the corresponding sites;

III. For AF=0.005%, the high-frequency sites (e.g. 5% of BRAF p.V600E) were used, and the frequency changes of the highfrequency sites were tested by ddPCR (e.g. BRAF p.V600E) diluted 100 fold is 0.05%) to determine the accuracy of the dilution fold, thereby *indirectly verifying* the frequency changes of the corresponding sites (for 0.05% of BRAF p.V600E was tested multiple times using optimized ddPCR, with the following requirements: the LOB of the site was controlled below 5 copies, the effective copy number was greater than 12,000, and the average value was calculated

IV. For AF=0%, the ddPCR conditions consistent with 0.5% were used, and the positive value was lower than the LOB value while ensuring that the number of copies was greater than

## **CONTACT US**